



HI-LO SECTIONAL TREATMENT TABLES

A. MODEL 6071



Common Features:

- “Hands-free” foot control operation.
- Open base allows closer access to patient.
- Powder coated heavy duty steel frame
- Quiet and stable, single column actuator adjusts from 24” to 33” H.
- Single column actuator saves space, no top shift when elevating.
- Hospital grade power cord.
- 100V-240V, 50/60Hz, 3.15 Amps



**A. Model 6071
2-Section Hi-Lo Treatment Table with
Armrests and Adjustable Glides**

- 2-Section Power Table can be used for a wide variety of purposes.
- Folding armrests for extra patient comfort.
- Includes Paper Dispenser (factory installed).
- 2-Section Top is 75” L x 28” W, with 2” thick high density foam.
- Head/Backrest is 25” L x 28” W & elevates up to 67° with dual release friction lock.
- Adjustable Glides allow for easy leveling.
- Agate Blue vinyl upholstery.

MODEL	LENGTH	WIDTH	HEIGHT
6071	75”	28”	24” - 33”

B. MODEL 6074



**B. Model 6074
3-Section, Mobile, Hi-Lo Treatment Table
with Locking Casters**

- 3-Section Power Table can be used for a wide variety of purposes.
- (4) 3” diameter, swiveling and locking casters.
- 3-Section Top is 78” L x 28” W with 2” thick high density foam.
- Contoured face/nose cutout in head section.
- Head section, 19” L x 28” W adjusts from flat to 90° below horizontal.
- Middle section, 21” L x 28” W is stationary.
- Leg/Backrest section 38” L x 28” W and elevates up to 70° with dual release friction lock.

MODEL	LENGTH	WIDTH	HEIGHT
6074	78”	28”	24” - 33”